

# LCGreen® Plus+

Gene Scanning Reagent

Information Sheet

*“The only dye for high-throughput, Hi-Res Melting™ analysis.”*



Our new LCGreen PLUS is manufactured exclusively by Idaho Technology and is designed specifically for high resolution melting curve analysis to detect DNA sequence variants (SNP's, insertions / deletions).

## Product Information

- LCGreen dyes are specifically designed for high-resolution melting curve analysis to detect DNA sequence variants (mutations, polymorphisms, etc).
- LCGreen PLUS is a new member of the dye family tailored for use in melting instruments with 96- or 384-well microtiter plates.
- LCGreen PLUS has superb fluorescence intensity, and can be used with other fluorescence based PCR detection systems such as the Roche LightCycler®. For optimal performance, the use of a high-resolution melting instrument is required.
- Optimum excitation: 440 – 470 nm. Optimum emission: 470 – 520 nm. Spectral characteristics depend on buffer composition, pH, ionic strength, and nucleic acid content of the solution.
- Addition of LCGreen dyes increase the melting temperature (T<sub>m</sub>) of DNA by about 1 – 3 °C, and may require adjustment of cycling parameters.
- LCGreen dyes are manufactured exclusively by Idaho Technology, and their chemical structures are unique among the scientific and patent literature. Patent pending.

## How it is Used

- LCGreen PLUS dye is supplied as a 10X solution in 10mM Tris-HCl, pH 7.4, 0.1 mM EDTA.
- LCGreen PLUS should be **used at 1X** for PCR. Add one volume of 10X solution to nine volumes of the PCR mixture.
- If you are using glass capillary tubes for PCR and/or for melting analysis, make sure your reaction mixture contains bovine serum albumin (BSA) at 250 - 500 µg/mL. BSA helps avoid enzyme, DNA and dye adhesion to the glass surface.


## Shipping & Storage

- Product is shipped at ambient temperature.
- Store at –20 °C upon receipt. Store at 4 °C after first use.
- Product is stable for one year at –20 °C, and up to 2 months at 4 °C.

## Related Products

- LightScanner™ (96 or 384 samples, high-resolution melting instrument)
- HR-1™ (single sample, high-resolution melting instrument)
- LCGreen® I gene scanning reagents (for use with Idaho Technology's HR-1)
- 10X BSA

**Conventional dsDNA dyes cannot be used at high concentrations due to dye redistribution during melting curve analysis.**



**Saturation of dsDNA binding sites eliminates potential for dye redistribution during melting curve acquisition.**

## Package Sizes

<b>No. of Reactions*</b>	1,000	10,000	Larger sizes Inquire
<b>LCGreen PLUS (10X solution)</b>	1 mL	10 X 1 mL	
<b>Catalog No.</b>	BCHM-ASY-0005	BCHM-ASY-0006	

\* based on 10 µl reaction volume


LCGreen, HR-1, LightScanner are trademarks of Idaho Technology Inc. LightCycler is a trademark of a member of the Roche Group. SYBR is a trademark of Molecular Probes

## References

1. Wittwer CT, GH Reed, CN Gundry, JG Vandersteen, and RJ Pryor. **High-Resolution Genotyping by Amplicon Melting Analysis using LC Green.** Clin. Chem., 49:853-60, 2003.
2. Liew M, Pryor R, Palais R, Meadows C, Erali M, Lyon E, Wittwer CT. **Genotyping of single nucleotide polymorphisms by high-resolution melting of small amplicons.** Clin. Chem., 50:1156-64, 2004.
3. Reed GH, Wittwer CT. **Sensitivity and specificity of SNP scanning by high-resolution melting analysis.** Clin. Chem., 50:1748-54, 2004.
4. Zhou L, Myers AN, Vandersteen JG, Wang L, Wittwer CT. **Closed-Tube Genotyping with Unlabeled Oligonucleotide Probes and a Saturating DNA Dye.** Clin Chem., 50:1328-35, 2004
5. McKinney JT, Longo N, Hahn SH, Matern D, Rinaldo P, Strauss AW, Dobrowolski SF. **Rapid, comprehensive screening of the human medium chain acyl-CoA dehydrogenase gene.** Mol Genet Metab. 82:112-20, 2004.
6. Willmore BS, Holden JA, Zhou L, Tripp S, Wittwer CT, Layfield LJ. **Detection of c-kit activating mutations in gastrointestinal stromal tumors by high-resolution amplicon melting analysis.** Am. J. Clin. Path., 122:206-16, 2004.
7. Zhou L Vandersteen J, Wang L, Fuller T, Taylor M, Palais B, Wittwer CT. **High-resolution DNA melting curve analysis to establish HLA genotypic identity.** Tissue Antigens, 64:156-64, 2004.
8. Vaughn CP, Elenitoba-Johnson KS. **High-resolution melting analysis for detection of internal tandem duplications.** J. Mol. Diag. 6:211-6, 2004.
9. Odell ID, Cloud JL Seipp M, and Wittwer CT. **Rapid species identification within the Mycobacterium chelonae-abscessus group by high-resolution melting of hsp65 PCR products.** Am J Clin Path 2005;123:96-101.
10. Dobrowolski SF, McKinney JT, Amat di San Filippo C, Giak Sim K, Wilcken B, Longo N. **Validation of dye-binding/high-resolution thermal denaturation for the identification of mutations in the SLC22A5 gene.** Hum Mutat. 2005 Mar;25(3):306-13

## Reviews on High-Resolution Melting

1. Gingeras TR, Higuchi R, Kricka LJ, Lo YM, Wittwer CT. **Fifty Years of Molecular (DNA/RNA) Diagnostics.** Clin Chem., 51:661-71, 2005
2. Highsmith WE Jr., **SNPs for sale. Cheap!.** Clin Chem., 50:1296-8, 2004



**DNA visualized with LCGreen PLUS dye in a 96-well plate with Idaho Technology's LightScanner™ instrument.**