

CAS-1200™

automated PCR setup

a Compact Robotic Workstation
for Precision Liquid Handling

www.PCRsetup.com

 **corbett**
LIFE SCIENCE

CAS-1200™

automated PCR setup


- ⊕ Addresses the precision demands of quantitative PCR setup
- ⊕ Compact size
- ⊕ Easy to use
- ⊕ Extremely versatile
- ⊕ Inexpensive to own and operate



← 40 cm (19") →

Why Automate PCR Set-Up?

- ✦ **Better reproducibility**
day-to-day, laboratory-to-laboratory
- ✦ **Reduce artefacts from human error**
- ✦ **Reduce Costs**
Use smaller reaction volumes
Fewer mistakes mean fewer repeat experiments
- ✦ **Reduced risk of repetitive strain injuries for staff**
- ✦ **Set-up is the least fun, but most critical part!**
The more experienced and careful you are the better the real-time PCR results—leaving it to a robot makes sense




3

CAS-1200™
automated PCR setup

EXAMPLE DATA

Real-time PCR analysis of robot-pipetted replicates

Lowest C_T 20.48
Highest C_T 20.66
 C_T Std Dev 0.05




β -actin samples from 12 wells set up in a 96-well plate (3 μ L sample with 17 μ L master mix)

4


Setup of a 2-fold Quantitative Amplification Series

- ✦ Human genomic DNA template diluted into 2-fold standards by the CAS-1200
- ✦ Amplification reactions were set up in replicates of four by the CAS-1200 using human dilutions as template
- ✦ Target: human BCL-2 gene
- ✦ Rotor-Gene 6000 real-time rotary analyzer was used for amplification and analysis
- ✦ Plots show tight replicates with expected one-cycle difference between each 2-fold dilution
- ✦ Standard curve for this experiment is shown with statistical data (boxed)




Comparison of manual vs. hand pipetting

20 µL reaction volume, 18 replicates



Hand pipetting

C_T std dev 0.07




CAS-1200 robot


C_T std dev 0.05

Comparison of manual vs. hand pipetting

15 μ L reaction volume, 18 replicates




Hand pipetting
C_T std dev 0.11




CAS-1200 robot
C_T std dev 0.07

Comparison of manual vs. hand pipetting

10 μ L reaction volume, 18 replicates



Hand pipetting
C_T std dev 0.12




CAS-1200 robot
C_T std dev 0.10

CAS-1200™ EXAMPLE DATA


automated PCR setup

Comparison of manual vs. hand pipetting

5 µL reaction volume, 18 replicates



Hand pipetting
C_T std dev 0.64



CAS-1200 robot
C_T std dev 0.12

CAS-1200™


automated PCR setup

Supports all amplification and real-time PCR formats


96-well and 384-well
(Incl. all ABI, Roche, Bio-Rad,
Stratagene, Eppendorf
etc. instruments)



SmartCycler™ Tubes



Rotor-Gene™ 0.2 mL, 0.1 mL
and Gene-Disc™ formats



LightCycler™ Capillaries



Roche COBAS™

CAS-1200™

automated PCR setup

96 and 384-well support




Compatible with plates, tubes, tube racks, cooling blocks etc

CAS-1200™

automated PCR setup

Unusual formats are also supported



Including: MALDI plates, ELISA plates, reservoir plates, etc, etc


CAS-1200™

automated PCR setup

www.corbettlifescience.com




Support for Rotor-Gene rotary tube formats




0.2 mL tubes



0.1 mL tubes



Gene-Disc™ 72



Gene-Disc™ 100

13

CAS-1200™

automated PCR setup

www.corbettlifescience.com



Direct setup of Rotor-Gene Gene-Disc™ tube format



Vertical tube orientation means a robot can set up all reactions directly into a Gene-Disc plate (both 72-well and 100-well discs)



Gene-Disc rotor ready for cycling

Heat Sealer provides a permanent or removable seal (user selectable)

14

Automatic normalization of sample concentrations

1. Specify the input concs. of each starting sample in the sample sheet
2. Specify the final desired conc.

3. Specify the final volume

Done!

(required vols. are specified)

Well	Contents
A1	118.7uL Water, 81.3uL Sample
B1	104.76uL Water, 95.24uL Sample B1 @ B2
C1	94.74uL Water, 105.26uL Sample C1 @ B2
D1	131.03uL Water, 68.97uL Sample D1 @ B2
E1	150.25uL Water, 49.75uL Sample A2 @ B2
F1	125.37uL Water, 74.63uL Sample B2 @ B2

Well	Sample Name	Sample ID	Conc.
A1	Sample A1 @ B2	Patient 1	123
B1	Sample B1 @ B2	Patient 2	105
C1	Sample C1 @ B2	Patient 3	95
D1	Sample D1 @ B2	Patient 4	145
A2	Sample A2 @ B2	Patient 5	201
B2	Sample B2 @ B2	Patient 6	134
C2	Sample C2 @ B2	Patient 7	176
D2	Sample D2 @ B2	Patient 8	190
A3	Sample A3 @ B2	Patient 9	98
B3	Sample B3 @ B2	Patient 10	143


Final Volume: Final Concentration:

Final Volume:

Final Volume:

Features

- ✦ Integrated cover
- ✦ Single channel pipetting
- ✦ UV light sterilisation option
- ✦ HEPA filter option
- ✦ Variety of plate configurations
- ✦ Variety of mixing possibilities
- ✦ Software wizards to setup standard and real-time PCR




CAS-1200™

automated PCR setup

Pipetting

- ✦ Single pipetting head
- ✦ Pipette volumes as low as 1 μL
- ✦ < 1% C.V. at 5 μL vol.
- ✦ Uses filtered Tecan or Capp style 50 μL and 200 μL tips
- ✦ Conductive tips are supported with optional liquid level sensing



Ready-racked Capp style tips




Tecan Genesis style conductive tips

CAS-1200™

automated PCR setup

Liquid level sensing (LLS)

- ✦ Optional circuitry
- ✦ Uses conductive tips (black, carbon impregnated)
- ✦ Liquid sensing as low as 10 μL
- ✦ Curbs tip immersion artifacts




CAS-1200™

automated PCR setup

Reagents and Master Mixes

- ✦ Removable Master Mix and Reagent block
- ✦ Can be refrigerated to keep reagents cool
- ✦ Many configurations available
- ✦ Custom versions available for master mix and reagent positions



19

Tip eject chute (to a bag or box)

1 × 200 μ L tip rack

Diluent (PCR water)

Master Mix Tubes

3 × 50 μ L tip racks

PCR Standards (set up by robot)

Reagents

DNA samples (input plate)

PCR reactions (output plate)

Typical deck layout for PCR set-up

CAS-1200™
automated PCR setup


20

CAS-1200™

automated PCR setup

Software

- ✦ Intuitive layout—simple to learn
- ✦ Standard PCR components included
- ✦ Prepare Master Mixes and controls
- ✦ Set-up of standards (dilution series)
- ✦ Normalize sample concentrations
- ✦ Tip re-use options
- ✦ Import/Export of sample names
- ✦ Pre- and Post run reports
- ✦ Virtual mode aids new protocol testing
- ✦ Most versatile software available



Robotic Workflow for Real-Time Analysis



X-tractor Gene™
automated nucleic acid extraction



CAS-1200™
automated sample setup



Rotor-Gene™
real-time rotary analyzer

“Extraction-to-Reaction”

Standard objections to robotics

- ⊕ **Luxury—only for big labs**
 - ◆ CAS-1200 costs less than a real-time instrument
 - ◆ Introduces a new generation of *personal* robotic workstation
- ⊕ **Running costs are too high**
 - ◆ Robotic tips are now similar in cost to regular tips
 - ◆ In fact, overall **savings** can be made on reagents & fewer repeat expts.
- ⊕ **Robots are too big—“we have no space”**
 - ◆ Only 49 cm (19.3”) wide—a fraction the size of similar equipment
- ⊕ **Difficult to use/can't be bothered learning to use it**
 - ◆ Try the software and see! We are happy to arrange a demonstration.

23


CAS-1200™ automated PCR setup

Summary

- ⊕ Low cost, small, simple to operate personal robot
- ⊕ High precision for real-time PCR setup
- ⊕ Low volume reaction setup saves reagent cost
- ⊕ Useful for many laboratory projects
- ⊕ Reduce human error and repetitive strain injury issues
- ⊕ Versatile; set up in 96-well, 384-well or wide range of other tube types
- ⊕ Single head for precision and cherry-picking of any sample
- ⊕ Optional Liquid Level Sensing, UV light sterilization & HEPA filter options



24



Offices

Brisbane Australia

Corbett Robotics Pty Ltd
42 McKechnie Drive
Eight Mile Plains, QLD 4113
T +61 7 3841 7077
F +61 7 3841 6077

Sydney Australia

Corbett Research Pty Ltd
14 Hilly Street
Mortlake, NSW 2137
T 1 800 803 915 (Toll free)
T +61 2 9736 1320
F +61 2 9736 1364

E-mail info@corbettlifescience.com

United Kingdom

Corbett Research UK Limited
Unit 1 Terek House, Sandpiper Court
Phoenix Business Park, Eaton, Socon
St Neots, Cambridgeshire, PE19 8EP
T +44 (0)1480 407 222
F +44 (0)1480 407 999

USA

Corbett Robotics Inc
185 Berry Street, Suite 5200
San Francisco, CA 94107
T +1 866 380 1166 (Toll free)
T +1 415 348 1166
F +1 415 348 1177

Web www.corbettlifescience.com

